

For Research Use Only

SLC25A28 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83127-5-PBS

Basic Information

Catalog Number:

83127-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34507

GenBank Accession Number:

BC076399

GeneID (NCBI):

81894

UNIPROT ID:

Q96A46

Full Name:

solute carrier family 25, member 28

Purification Method:

Protein A purification

CloneNo.:

230392F8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00142-2, SLC25A28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83127-2-PBS. Detection antibody: 83127-5-PBS. Standard: Ag34507. Range: 1.25-20 ng/mL.