

For Research Use Only

GRIA2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83128-2-PBS



Basic Information

Catalog Number:

83128-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2662

GenBank Accession Number:

BC010574

GeneID (NCBI):

2891

UNIPROT ID:

P42262

Full Name:

glutamate receptor, ionotropic, AMPA

2

Calculated MW:

883 aa, 99 kDa

Purification Method:

Protein A purification

CloneNo.:

230403D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00143-1, Glutamate receptor 2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83128-2-PBS. Detection antibody: 83128-1-PBS. Standard: Ag2662. Range: 1.25-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83128-2-PBS against Human GRIA2 was performed. The affinity constant is 5.96 nM.