

For Research Use Only

NLRP10 Recombinant antibody, PBS Only (Capture)

Catalog Number:83131-1-PBS



Basic Information

Catalog Number: 83131-1-PBS	GenBank Accession Number: BC104957	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 338322	CloneNo.: 230402C6
Source: Rabbit	UNIPROT ID: Q86W26	
Isotype: IgG	Full Name: NLR family, pyrin domain containing 10	
Immunogen Catalog Number: AG26402	Calculated MW: 655 aa, 75 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00242-2, NLRP10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83131-1-PBS. Detection antibody: 83131-7-PBS. Standard: Ag26402. Range: 2.5-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83131-1-PBS against Human NLRP10 was performed. The affinity constant is below 1 pM.