

For Research Use Only

MTCO2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83132-4-PBS

Basic Information

Catalog Number:

83132-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

YP_003024029

GeneID (NCBI):

4513

Full Name:

cytochrome c oxidase II

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

230448C11

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00254-3, MTCO2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83132-4-PBS. Detection antibody: 83132-6-PBS. Standard: SY01099. Range: 0.313-20 ng/mL