

For Research Use Only

# SPTA1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83133-7-PBS



## Basic Information

Catalog Number:

83133-7-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_003126

GeneID (NCBI):

6708

UNIPROT ID:

P02549

Full Name:

spectrin, alpha, erythrocytic 1  
(elliptocytosis 2)

Calculated MW:

280 kDa

Purification Method:

Protein A purification

CloneNo.:

230507G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00280-2, SPTA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83133-7-PBS. Detection antibody: 83133-5-PBS. Standard: SY01287. Range: 78-2500 pg/mL.