

For Research Use Only

# JUND Recombinant antibody, PBS Only

Catalog Number: 83134-4-PBS



## Basic Information

<b>Catalog Number:</b> 83134-4-PBS	<b>GenBank Accession Number:</b> NM_005354	<b>Purification Method:</b> Protein A purification
<b>Concentration:</b> 1mg/ml	<b>GeneID (NCBI):</b> 3727	<b>CloneNo.:</b> 230457E8
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P17535	
<b>Isotype:</b> IgG	<b>Full Name:</b> jun D proto-oncogene	
	<b>Calculated MW:</b> 35 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA  
**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data