

For Research Use Only

ORM2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83135-1-PBS



Basic Information

Catalog Number:

83135-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1667

GenBank Accession Number:

BC015964

GeneID (NCBI):

5005

UNIPROT ID:

P19652

Full Name:

orosomuroid 2

Calculated MW:

24 kDa

Purification Method:

Protein A purification

CloneNo.:

230477B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00283-2, ORM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83135-1-PBS. Detection antibody: 83135-3-PBS. Standard: Ag1667. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00283-3, ORM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83135-1-PBS. Detection antibody: 83135-2-PBS. Standard: Ag1667. Range: 1.25-20 ng/mL



Bilayer interferometry (BLI) kinetic assay of 83135-1-PBS against Human ORM2 was performed. The affinity constant is 0.114 nM.