

For Research Use Only

RNF31 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83137-1-PBS

Basic Information

Catalog Number: 83137-1-PBS	GenBank Accession Number: BC017376	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 55072	CloneNo.: 230446B3
Source: Rabbit	UNIPROT ID: Q96EPO	Affinity: $K_D = 2.69 \times 10^{-9}$
Isotype: IgG	Full Name: ring finger protein 31	$K_{off} = 3.63 \times 10^{-4}$
Immunogen Catalog Number: AG34505	Calculated MW: 1072 aa, 120 kDa	$K_{on} = 1.35 \times 10^5$

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00285-1, RNF31 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83137-3-PBS. Detection antibody: 83137-1-PBS. Standard: Ag34505. Range: 1.25-20 ng/mL.