

For Research Use Only

# CTNND2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83138-2-PBS

## Basic Information

Catalog Number:

83138-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001332

GeneID (NCBI):

1501

UNIPROT ID:

Q9UQB3

Full Name:

catenin (cadherin-associated protein),  
delta 2 (neural plakophilin-related  
arm-repeat protein)

Calculated MW:

133 kDa

Purification Method:

Protein A purification

CloneNo.:

230396D12

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00295-2, CTNND2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83138-1-PBS. Detection antibody: 83138-2-PBS. Standard: SY00922. Range: 0.078-2.5 ng/mL.