

For Research Use Only

CTNND2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83138-3-PBS

Basic Information

Catalog Number:

83138-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001332

GeneID (NCBI):

1501

UNIPROT ID:

Q9UQB3

Full Name:

catenin (cadherin-associated protein),
delta 2 (neural plakophilin-related
arm-repeat protein)

Calculated MW:

133 kDa

Purification Method:

Protein A purification

CloneNo.:

230396D7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00295-1, CTNND2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83138-4-PBS. Detection antibody: 83138-3-PBS. Standard: SY00922. Range: 0.078-2.5 ng/mL