

For Research Use Only

PPFIA3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83142-1-PBS

Basic Information

Catalog Number:

83142-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1971

GenBank Accession Number:

BC021255

GeneID (NCBI):

8541

UNIPROT ID:

O75145

Full Name:

protein tyrosine phosphatase,
receptor type, f polypeptide (PTPRF),
interacting protein (liprin), alpha 3

Calculated MW:

1194 aa, 133 kDa

Purification Method:

Protein A purification

CloneNo.:

230459A9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00170-1, PPFIA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83142-1-PBS. Detection antibody: 83142-3-PBS. Standard: Ag1971. Range: 1.25-20 ng/mL.