

For Research Use Only

# TMEM165 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83146-2-PBS

## Basic Information

Catalog Number:

83146-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14362

GenBank Accession Number:

BC003545

GeneID (NCBI):

55858

UNIPROT ID:

Q9HC07

Full Name:

transmembrane protein 165

Calculated MW:

324 aa, 35 kDa

Purification Method:

Protein A purification

CloneNo.:

230517C10

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve ofMP00217-1, TMEM165 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83146-2-PBS. Detection antibody: 83146-3-PBS. Standard: Ag14362. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00217-2, TMEM165 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83146-2-PBS. Detection antibody: 83146-1-PBS. Standard: Ag14362. Range: 0.625-20 ng/mL