

For Research Use Only

GRP94 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83147-3-PBS

Basic Information

Catalog Number: 83147-3-PBS	GenBank Accession Number: BC009195	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 7184	CloneNo.: 230502F1
Source: Rabbit	UNIPROT ID: P14625	Affinity: $K_D < 1.00 \times 10^{-12}$
Isotype: IgG	Full Name: heat shock protein 90kDa beta (Grp94), member 1	$K_{off} < 1.00 \times 10^{-6}$ $K_{on} = 1.53 \times 10^5$
Immunogen Catalog Number: AG1439	Calculated MW: 96 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00194-2, GRP94 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83147-2-PBS. Detection antibody: 83147-3-PBS. Standard: Ag1439. Range: 0.625-20 ng/mL