

For Research Use Only

# IL-17RB Recombinant antibody, PBS Only (Detector)

Catalog Number: 83149-1-PBS



## Basic Information

Catalog Number:

83149-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_018725

GeneID (NCBI):

55540

UNIPROT ID:

Q9NRM6

Full Name:

interleukin 17 receptor B

Calculated MW:

56 kDa

Purification Method:

Protein A purification

CloneNo.:

230439E1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00240-1, IL-17RB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83149-2-PBS. Detection antibody: 83149-1-PBS. Standard: SY00888. Range: 0.078-2.5 ng/mL