

For Research Use Only

# IL-17RB Recombinant antibody

Catalog Number: 83149-3-RR



## Basic Information

Catalog Number:

83149-3-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_018725

GeneID (NCBI):

55540

UNIPROT ID:

Q9NRM6

Full Name:

interleukin 17 receptor B

Calculated MW:

56 kDa

Purification Method:

Protein A purification

CloneNo.:

230439B3

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

## Background Information

IL17RB, also named as CRL4, EVI27 and IL17BR, is a receptor for the proinflammatory cytokines IL17B and IL17E. IL17RB may play a role in controlling the growth and/or differentiation of hematopoietic cells. IL17RB has two isoforms with MW 56 kDa (isoform 1) and 32 kDa (isoform 2). The antibody is specific to IL17RB.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

