

For Research Use Only

SLC6A9 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83150-3-PBS



Basic Information

Catalog Number:

83150-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34437

GenBank Accession Number:

NM_201649

GeneID (NCBI):

6536

UNIPROT ID:

P48067

Full Name:

solute carrier family 6
(neurotransmitter transporter,
glycine), member 9

Calculated MW:

706 aa, 78 kDa

Purification Method:

Protein A purification

CloneNo.:

240010H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00243-1, SLC6A9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83150-3-PBS. Detection antibody: 83150-1-PBS. Standard: Ag34437. Range: 0.625-20 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83150-3-PBS against Human SLC6A9 was performed. The affinity constant is 0.144 nM.