

For Research Use Only

TIMP1 Recombinant antibody, PBS Only

Catalog Number: 83151-6-PBS



Basic Information

Catalog Number:

83151-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC000866

GeneID (NCBI):

7076

UNIPROT ID:

P01033

Full Name:

TIMP metalloproteinase inhibitor 1

Calculated MW:

23 kDa

Purification Method:

Protein A purification

CloneNo.:

230329D8

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

