

For Research Use Only

B7-H3/CD276 Recombinant antibody, PBS Only

Catalog Number: 83152-5-PBS



Basic Information

Catalog Number:

83152-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001024736

GeneID (NCBI):

80381

UNIPROT ID:

Q5ZPR3

Full Name:

CD276 molecule

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

230305A10

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83152-5-RR (B7-H3/CD276 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83152-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83152-5-RR against Human CD276 were performed. The affinity constant is 0.607 nM.