

For Research Use Only

FMO1 Recombinant antibody, PBS Only



Catalog Number: 83156-3-PBS

Basic Information

Catalog Number:

83156-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34267

GenBank Accession Number:

NM_001282692

GeneID (NCBI):

2326

UNIPROT ID:

Q01740

Full Name:

flavin containing monooxygenase 1

Calculated MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

230455B11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83156-3-RR (FMO1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83156-3-PBS in a different storage buffer formulation.