

For Research Use Only

Mouse MCP-1/CCL2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83157-2-PBS



Basic Information

Catalog Number:

83157-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM-011333

GeneID (NCBI):

20296

Full Name:

chemokine (C-C motif) ligand 2

Purification Method:

Protein A purification

CloneNo.:

230483G3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00168-1, Mouse Ccl2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83157-4-PBS. Detection antibody: 83157-2-PBS. Standard: Eg0427. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00168-2, Mouse Ccl2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83157-3-PBS. Detection antibody: 83157-2-PBS. Standard: Eg0427. Range: 0.625-20 ng/mL