

For Research Use Only

IFI35 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83160-3-PBS



Basic Information

Catalog Number:

83160-3-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34653

GenBank Accession Number:

BC001356

GeneID (NCBI):

3430

UNIPROT ID:

P80217

Full Name:

interferon-induced protein 35

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

240063H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00191-1, IFI35 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83160-3-PBS. Detection antibody: 83160-1-PBS. Standard: Ag34653. Range: 0.625-20 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83160-3-PBS against Human IFI35 was performed. The affinity constant is 0.343 nM.