

For Research Use Only



TYK2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83161-3-PBS

Basic Information

Catalog Number: 83161-3-PBS	GenBank Accession Number: BC014243	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 7297	CloneNo.: 230176F8
Source: Rabbit	UNIPROT ID: P29597	
Isotype: IgG	Full Name: tyrosine kinase 2	
Immunogen Catalog Number: AG28457	Calculated MW: 1187 aa, 134 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00203-2, TYK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83161-1-PBS. Detection antibody: 83161-3-PBS. Standard: Ag28457. Range: 1.25-20 ng/mL.