

For Research Use Only

# HIC2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83162-1-PBS

## Basic Information

Catalog Number: 83162-1-PBS	GenBank Accession Number: BC094787	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 23119	CloneNo.: 230177A3
Source: Rabbit	UNIPROT ID: Q96JB3	
Isotype: IgG	Full Name: hypermethylated in cancer 2	
Immunogen Catalog Number: AG18789	Calculated MW: 615 aa, 66 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00206-1, HIC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83162-1-PBS. Detection antibody: 83162-5-PBS. Standard: Ag18789. Range: 0.625-20 ng/mL