

For Research Use Only

TSG-6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83163-2-PBS



Basic Information

Catalog Number:

83163-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3885

GenBank Accession Number:

BC030205

GeneID (NCBI):

7130

UNIPROT ID:

P98066

Full Name:

tumor necrosis factor, alpha-induced protein 6

Calculated MW:

277 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

230438E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00218-1, TSG-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83163-2-PBS. Detection antibody: 83163-3-PBS. Standard: Ag3885. Range: 0.625-20 ng/mL