

For Research Use Only

COX6B2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83164-2-PBS

Basic Information

Catalog Number: 83164-2-PBS	GenBank Accession Number: BC064548	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 125965	CloneNo.: 230497D6
Source: Rabbit	UNIPROT ID: Q6YFQ2	
Isotype: IgG	Full Name: cytochrome c oxidase subunit VIb polypeptide 2 (testis)	
Immunogen Catalog Number: AG1983	Calculated MW: 88 aa, 11 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00220-1, COX6B2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83164-2-PBS. Detection antibody: 83164-1-PBS. Standard: Ag1983. Range: 2.5-20 ng/mL