

For Research Use Only

# SAR1A Recombinant antibody, PBS Only (Capture)



Catalog Number: 83168-5-PBS

## Basic Information

Catalog Number:

83168-5-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_020150

GeneID (NCBI):

56681

UNIPROT ID:

Q9NR31

Full Name:

SAR1 homolog A (*S. cerevisiae*)

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

230471D4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00263-3, SAR1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83168-5-PBS. Detection antibody: 83168-4-PBS. Standard: SY01131. Range: 0.156-5 ng/mL.