

For Research Use Only

PGF/PIGF Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83170-2-PBS



Basic Information

Catalog Number:

83170-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0969

GenBank Accession Number:

BC007789

GeneID (NCBI):

5228

UNIPROT ID:

P49763

Full Name:

placental growth factor

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

230532C12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00146-1, PGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83170-4-PBS. Detection antibody: 83170-2-PBS. Standard: Ag0969. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00146-2, PGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83170-3-PBS. Detection antibody: 83170-2-PBS. Standard: Ag0969. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00146-3, PGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83170-2-PBS. Detection antibody: 83170-1-PBS. Standard: Ag0969. Range: 0.625-20 ng/mL