

For Research Use Only

# PGF/PIGF Recombinant antibody, PBS Only (Capture)

Catalog Number: 83170-3-PBS



## Basic Information

Catalog Number:

83170-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0969

GenBank Accession Number:

BC007789

GeneID (NCBI):

5228

UNIPROT ID:

P49763

Full Name:

placental growth factor

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

230532A9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00146-2, PGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83170-3-PBS. Detection antibody: 83170-2-PBS. Standard: Ag0969. Range: 0.625-20 ng/mL