

For Research Use Only

Histone H2B Recombinant antibody, PBS Only (Detector)



Catalog Number: 83171-3-PBS

Basic Information

Catalog Number:

83171-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC005827

GeneID (NCBI):

8349

UNIPROT ID:

Q16778

Full Name:

histone cluster 2, H2be

Calculated MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

230539H2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00183-1, Histone H2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83171-2-PBS. Detection antibody: 83171-3-PBS. Standard: SY02145Ac. Range: 0.625-5 ng/mL



Cytometric bead array standard curve ofMP00183-2, Histone H2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83171-1-PBS. Detection antibody: 83171-3-PBS. Standard: SY02145Ac. Range: 0.625-5 ng/mL