

For Research Use Only

ZGLP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83176-3-PBS



Basic Information

Catalog Number:

83176-3-PBS

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001103167

GeneID (NCBI):

100125288

UNIPROT ID:

POC6A0

Full Name:

GATA like protein-1

Calculated MW:

30 kDa

Purification Method:

Protein A purification

CloneNo.:

230407E12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00250-1, ZGLP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83176-3-PBS. Detection antibody: 83176-2-PBS. Standard: SY01346. Range: 0.078-1.25 ng/mL