

For Research Use Only

Rat IL-6 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83179-1-PBS

Basic Information

Catalog Number:

83179-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_012589.2

GeneID (NCBI):

24498

UNIPROT ID:

P20607

Full Name:

interleukin 6

Calculated MW:

24 kDa

Purification Method:

Protein A purification

CloneNo.:

240003B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00264-1, Rat IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83179-1-PBS. Detection antibody: 83179-3-PBS. Standard: Eg0800. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00264-3, Rat IL-6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83179-2-PBS. Detection antibody: 83179-1-PBS. Standard: Eg0800. Range: 0.625-20 ng/mL