

For Research Use Only

MKS1 Recombinant antibody

Catalog Number: 83181-2-RR



Basic Information

Catalog Number:

83181-2-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9177

GenBank Accession Number:

BC010061

GeneID (NCBI):

54903

UNIPROT ID:

Q9NXB0

Full Name:

Meckel syndrome, type 1

Calculated MW:

559 aa, 65 kDa

Purification Method:

Protein A purification

CloneNo.:

240036E12

Applications

Tested Applications:

ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83181-2-RR against Human BBS13/MKS1 were performed. The affinity constant is 0.198 nM.