

For Research Use Only

KLF9 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83183-2-PBS

Basic Information

Catalog Number: 83183-2-PBS	GenBank Accession Number: BC074880	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 687	CloneNo.: 230436G4
Source: Rabbit	UNIPROT ID: Q13886	
Isotype: IgG	Full Name: Kruppel-like factor 9	
Immunogen Catalog Number: AG34411	Calculated MW: 27 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00148-2, KLF9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83183-2-PBS. Detection antibody: 83183-4-PBS. Standard: Ag34411. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00148-3, KLF9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83183-2-PBS. Detection antibody: 83183-1-PBS. Standard: Ag34411. Range: 1.25-20 ng/mL