

For Research Use Only

JAM-A/CD321 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83186-3-PBS



Basic Information

Catalog Number:

83186-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC001533

GeneID (NCBI):

50848

UNIPROT ID:

Q9Y624

Full Name:

F11 receptor

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

240051H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00247-1, F11R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83186-2-PBS. Detection antibody: 83186-3-PBS. Standard: SY00088. Range: 0.625-20 ng/mL