

For Research Use Only

SMC6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83187-4-PBS

Basic Information

Catalog Number:

83187-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6898

GenBank Accession Number:

BC039828

GeneID (NCBI):

79677

UNIPROT ID:

Q96SB8

Full Name:

structural maintenance of chromosomes 6

Calculated MW:

126 kDa

Purification Method:

Protein A purification

CloneNo.:

240015E9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00255-3, SMC6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83187-4-PBS. Detection antibody: 83187-5-PBS. Standard: Ag6898. Range: 1.25-20 ng/mL.