

For Research Use Only

# MMP13 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83188-1-PBS

## Basic Information

Catalog Number:

83188-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34892

GenBank Accession Number:

BC074808

GeneID (NCBI):

4322

UNIPROT ID:

P45452

Full Name:

matrix metalloproteinase 13

Calculated MW:

471 aa, 54 kDa

Purification Method:

Protein A purification

CloneNo.:

240060A6

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00266-2, Human MMP13  
Recombinant Matched Antibody Pair - BSA and  
Azide Free measured by Cytometric bead array.  
Capture antibody: 83188-1-PBS. Detection  
antibody: 83188-2-PBS. Standard: Ag34892. Range:  
5-40 ng/mL.