

For Research Use Only

MMP13 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83188-3-PBS

Basic Information

Catalog Number: 83188-3-PBS	GenBank Accession Number: BC074808	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4322	CloneNo.: 240060D8
Source: Rabbit	UNIPROT ID: P45452	
Isotype: IgG	Full Name: matrix metalloproteinase 13	
Immunogen Catalog Number: AG34892	Calculated MW: 471 aa, 54 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00266-1, MMP13 Recombinant Matched Antibody Pair, PBS Only. 83188-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag34892. 83188-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



Cytometric bead array standard curve of MP00266-1, MMP13 Recombinant Matched Antibody Pair, PBS Only. PBS Only. Capture antibody: 83188-3-PBS. Detection antibody: 83188-2-PBS. Standard: Ag34892. Range: 5-40 ng/mL