

For Research Use Only

FOXP2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83190-1-PBS

Basic Information

Catalog Number:	83190-1-PBS	GenBank Accession Number:	BC018016	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	93986	CloneNo.:	240093B2
Source:	Rabbit	UNIPROT ID:	O15409		
Isotype:	IgG	Full Name:	forkhead box P2		
Immunogen Catalog Number:	AG16559	Calculated MW:	80 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00273-1, FOXP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83190-2-PBS. Detection antibody: 83190-1-PBS. Standard: Ag16559. Range: 2.5-40 ng/mL.