

For Research Use Only

PROCA1 Recombinant antibody, PBS Only

Catalog Number: 83193-3-PBS



Basic Information

Catalog Number:

83193-3-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10851

GenBank Accession Number:

BC029574

GeneID (NCBI):

147011

UNIPROT ID:

Q8NCQ7

Full Name:

proline-rich cyclin A1-interacting protein

Calculated MW:

336 aa, 38 kDa

Purification Method:

Protein A purification

CloneNo.:

230538D11

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

