

For Research Use Only

# TFAP4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83196-2-PBS

## Basic Information

Catalog Number: 83196-2-PBS	GenBank Accession Number: BC010576	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 7023	CloneNo.: 240031G6
Source: Rabbit	UNIPROT ID: Q01664	
Isotype: IgG	Full Name: transcription factor AP-4 (activating enhancer binding protein 4)	
Immunogen Catalog Number: AG2633	Calculated MW: 338 aa, 39 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00288-1, TFAP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83196-3-PBS. Detection antibody: 83196-2-PBS. Standard: Ag2633. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00288-3, TFAP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83196-1-PBS. Detection antibody: 83196-2-PBS. Standard: Ag2633. Range: 0.625-20 ng/mL