

For Research Use Only

GPR35 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83197-1-PBS

Basic Information

Catalog Number: 83197-1-PBS	GenBank Accession Number: NM_005301	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2859	CloneNo.: 240032A11
Source: Rabbit	UNIPROT ID: Q9HC97	
Isotype: IgG	Full Name: G protein-coupled receptor 35	
	Calculated MW: 34 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00286-1, GPR35 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83197-1-PBS. Detection antibody: 83197-3-PBS. Standard: SY01275. Range: 78.125-5000 pg/mL