

For Research Use Only

GPR35 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83197-3-PBS

Basic Information

Catalog Number:

83197-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005301

GeneID (NCBI):

2859

UNIPROT ID:

Q9HC97

Full Name:

G protein-coupled receptor 35

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

240032B5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00286-1, GPR35 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83197-1-PBS. Detection antibody: 83197-3-PBS. Standard: SY01275. Range: 78.125-5000 pg/mL