

For Research Use Only

# GPR35 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83197-4-PBS



## Basic Information

Catalog Number:

83197-4-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005301

GeneID (NCBI):

2859

UNIPROT ID:

Q9HC97

Full Name:

G protein-coupled receptor 35

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

240032H4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00286-2, GPR35 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83197-2-PBS. Detection antibody: 83197-4-PBS. Standard: SY01275. Range: 78.125-5000 pg/mL