

For Research Use Only

Mouse MIF Recombinant antibody, PBS Only (Detector)

Catalog Number: 83199-1-PBS



Basic Information

Catalog Number:

83199-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_010798.2

GeneID (NCBI):

17319

UNIPROT ID:

P34884

Full Name:

macrophage migration inhibitory factor

Calculated MW:

13KD

Purification Method:

Protein A purification

CloneNo.:

230529C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00282-1, Mouse MIF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83199-2-PBS. Detection antibody: 83199-1-PBS. Standard: Eg0574. Range: 3.13-100 ng/mL