

For Research Use Only

# IGFBP4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83200-2-PBS

## Basic Information

Catalog Number:

83200-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001552.3

GeneID (NCBI):

3487

UNIPROT ID:

P22692-1

Full Name:

insulin-like growth factor binding protein 4

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

230543E12

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00279-1, Human IGFBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83200-2-PBS. Detection antibody: 83200-1-PBS. Standard: Eg0825. Range: 1.56-50 ng/mL