

For Research Use Only

Sestrin 2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83202-2-PBS



Basic Information

Catalog Number:

83202-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1247

GenBank Accession Number:

BC013304

GeneID (NCBI):

83667

UNIPROT ID:

P58004

Full Name:

sestrin 2

Calculated MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

240034G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00278-1, Sestrin 2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83202-2-PBS. Detection antibody: 83202-3-PBS. Standard: Ag1247. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00278-2, Sestrin 2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83202-2-PBS. Detection antibody: 83202-1-PBS. Standard: Ag1247. Range: 1.25-20 ng/mL