

For Research Use Only

CEACAM6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83203-3-PBS

Basic Information

Catalog Number:

83203-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC005008

GeneID (NCBI):

4680

UNIPROT ID:

P40199

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 6 (non-
specific cross reacting antigen)

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

230475C6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00144-2, CEACAM6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83203-3-PBS. Detection antibody: 83203-1-PBS. Standard: Eg0056. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00144-3, CEACAM6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83203-3-PBS. Detection antibody: 83203-2-PBS. Standard: Eg0056. Range: 0.625-10 ng/mL