

For Research Use Only

CEACAM6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83203-6-PBS



Basic Information

Catalog Number:

83203-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC005008

GeneID (NCBI):

4680

UNIPROT ID:

P40199

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 6 (non-
specific cross reacting antigen)

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

230475F2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP00144-4, Human CEACAM6 Monoclonal Matched Antibody Pair - PBS only. 83203-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0056. 83203-6-PBS was HRP conjugated as the detection antibody. Range: 7.8-500 pg/mL.



For the PBMC supernatant cultured for 3 days, the mean human CEACAM6 concentration was determined to be 513.9 pg/mL in PHA stimulated PBMC supernatant.



Serum of sixteen individual healthy human donors were measured. The human CEACAM6 concentration of detected samples was determined to be 11,882.1 pg/mL with a range of 4,674.3 - 21,548.7 pg/mL.