

For Research Use Only

KLF5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83205-3-PBS



Basic Information

Catalog Number:

83205-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15267

GenBank Accession Number:

BC042131

GeneID (NCBI):

688

UNIPROT ID:

Q13887

Full Name:

Kruppel-like factor 5 (intestinal)

Calculated MW:

457 aa, 51 kDa

Purification Method:

Protein A purification

CloneNo.:

230509D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00155-1, KLF5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83205-3-PBS. Detection antibody: 83205-1-PBS. Standard: Ag15267. Range: 0.625-20 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83205-3-PBS against Human KLF5 was performed. The affinity constant is 0.219 nM.