

For Research Use Only

# GRAF Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 83206-2-PBS



## Basic Information

Catalog Number:

83206-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12008

GenBank Accession Number:

BC068555

GeneID (NCBI):

23092

UNIPROT ID:

Q9UNA1

Full Name:

Rho GTPase activating protein 26

Calculated MW:

759 aa, 86 kDa

Purification Method:

Protein A purification

CloneNo.:

240086B4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00157-1, ARHGAP26 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83206-2-PBS. Detection antibody: 83206-1-PBS. Standard: Ag12008. Range: 0.625-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83206-2-PBS against Human GRAF was performed. The affinity constant is 66.5 nM.